

## **ACCESSWALK® Specifications & Warranty**

- Dimensions of ACCESSWALK® shipping in a cardboard box: 39" × 11" × 38" for 40 lb. (98 cm x 28 cm x 95.5 cm for 18 kg)
- ACCESSWALK® fits through standard ADA doorways
- Lightweight anodized aluminum frame
- Foldable for storage and transportation
- Folded Up Dimensions: (without 12" front wheels & back bar): 33.5" L (85.1 cm) x 34.45" W (87.5 cm) x 10" H (25.4 cm)
- Folded Up Dimensions: (with 12" front wheels): 33.5" L (85.1 cm) x 34.45" W (87.5 cm) x 17" H (43.2 cm)
- Width of ACCESSWALK®: 34.45" (87.5 cm)
- Length: 30" (76 cm)
- Net Weight: 37.5 lb. (17 kg)
- Gross Weight: 39.7 lb. (18 kg)
- Maximum weight capacity: 285 lb. (130 kg)
- Folding seat and textile storage basket
- Seat Material: Polyester mesh incl. foamed cushion
- Seat Height: 23" (59 cm)
- Seat Dimensions: 16.3" wide x 11.8" deep (41.5 cm x 30 cm)
- Parking brake: Yes, activated with the break hand brake
- Turning Radius Circle: 35.43" (90 cm)
- Front Wheel Stainless Steel casters/axel
- Yellow "SUNNY" Pneumatic rubber wheels: Two front wheels of 12" diameter (30.5 cm) and two back wheels of 10" diameter (25.4 cm)

### **ACCESSREC**

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: [sales@AccessRec.com](mailto:sales@AccessRec.com)

Web: [www.ACCESSRECC.com](http://www.ACCESSRECC.com)

- Adjustable/removable handles from 32.5" to 37" high (82.55 cm to 94 cm)
- Width Between Handles: 18" (45.7 cm)
- Handles material: Plastic
- Carry Basket: 11 lb. payload (5 kg)
- Carry Basket Material: Mesh fabric
- Warranty: Limited Warranty of 3 Years

**ACCESSREC**

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: [sales@AccessRec.com](mailto:sales@AccessRec.com)

Web: [www.ACCESSREcc.com](http://www.ACCESSREcc.com)



## **ACCESSREC Rollator Walker – Limited Warranty**

*ACCESSREC provides three (3) years limited warranty from the date of purchase, covering defects in materials and workmanship for the rollator frame only.*

### **Warranty Exclusions**

*The following components are not covered under this warranty:*

*Seat cover material*

*Carry basket*

*Brakes*

*Plastic handles or other plastic spare parts*

*Wheels*

*All warranty claims and limitations will be assessed solely by ACCESSREC.*

### **Warranty Limitations**

*The three (3) years limited warranty does not cover damage to the rollator frame caused by:*

*Normal wear and tear*

*Improper use, maintenance, or storage*

*Accidents or natural disasters*

*Unauthorized modifications made without prior written consent from ACCESSREC*

*Unauthorized repairs performed without prior written approval from ACCESSREC*

*For any warranty inquiries or claims, please contact ACCESSREC directly.*

### **ACCESSREC**

67 Sand Park Road - Suite A, Cedar Grove, New Jersey 07009 - USA

Ph: +1 973-955-0514 - Fax: +1 201-624-7007

Email: [sales@AccessRec.com](mailto:sales@AccessRec.com)

Web: [www.ACCESSRECC.com](http://www.ACCESSRECC.com)